

PlasmidExpress Nanopore Sequencing

\$4.99/sample

5 Free Samples per Account | Promo code: **5FREE**

Ideal alternative to standard whole plasmid sequencing, providing fast and convenient sequencing for plasmids up to 10kb. This service is available at a cost-effective rate of just \$4.99 per sample with overnight turnaround.



Industry-leading Technology

Nanopore sequencing combined with Quintara's proprietary technology



Affordable Price

\$4.99/sample flat rate for all sample types, no hidden fees



6 Types of Samples

- Plasmid DNA(10μL+, 50-100ng/μL)
- Bacterial Overnight Culture(10μL+)
- RCA Product(10μL+)
- Colony Suspension(10μL+)
- Agar Plate
- Glycerol Stock(10μL+)



Fast Turnaround

Overnight data delivery for all sample types



Easy submission

No primers needed




Every customer gets 5 free samples

For detailed information, including promotions, logistics, sample requirements and submission instructions, please go to our website

www.quintarabio.com/nanopore-sequencing-service

or email sales.us@quintarabio.com


Workflow



Sample Submission



Library Preparation



Nanopore Sequencing



Data Analysis

(FASTA File, Per-Base-Details File, Chromatogram File, Genbank Annotation File)

MASSACHUSETTS

E sales.us@quintarabio.com

T 1-617-943-2768

A 625 Mt Auburn St, Ste 105

CALIFORNIA

E sales.us@quintarabio.com

T 1-510-397-0215

A 3563 Investment Blvd, #2