

# Whole Plasmid Sequencing - Express

\$4.99/sample

Ideal alternative to standard whole plasmid sequencing, providing fast and convenient sequencing for plasmids up to 10kb. This service is available at a cost-effective rate of just \$4.99 per sample with overnight turnaround.



## Industry-leading Technology

Nanopore sequencing combined with Quintara's proprietary technology



## Affordable Price

\$4.99/sample flat rate for all sample types, no hidden fees



## 6 Types of Samples

- Plasmid DNA(10µL+, 50-100ng/µL)
- Bacterial Overnight Culture(10µL+)
- RCA Product(10µL+)
- Colony Suspension(10µL+)
- Agar Plate
- Glycerol Stock(10µL+)



## Fast Turnaround

Overnight data delivery for all sample types



## Easy submission

No primers needed




Receive a \$5 Amazon gift card when you register a Quintara account

For detailed information, including promotions, logistics, sample requirements and submission instructions, please go to our website

[www.quintarabio.com/nanopore-sequencing-service](http://www.quintarabio.com/nanopore-sequencing-service)

or email [sales.us@quintarabio.com](mailto:sales.us@quintarabio.com)


## Workflow



### Sample Submission



### Library Preparation



### Nanopore Sequencing



### Data Analysis

(FASTA File, Per-Base-Details File, Chromatogram File, Genbank Annotation File)

#### CA Hayward

T 1-510-397-0215

A 3563 Investment Blvd, #2  
Hayward, CA 94545

#### Massachusetts

T 1-617-943-2768

A 625 Mt Auburn Street, Suite 105  
Cambridge, MA 02138

#### Texas

T 346-239-9971

A 10410 Corporate Dr. Suite 118  
Sugar Land, TX 77478

#### Michigan

T 734-476-9826

A 600 S Wagner Rd, Suite 35  
Ann Arbor, MI 48103

#### CA San Diego

T 619-823-0895

A 6364 Ferris Square, Suite 102  
San Diego, CA 92121