

# Bispecific Antibody Service

Multiple types | Tailored services | Fast turnaround

## ▶ Quintara Accelerates Your Antibody Discovery

A bispecific antibody (BsAb) is a specialized type of antibody engineered to simultaneously bind two distinct antigens or different epitopes on the same antigen. It is a powerful tool in cancer immunotherapy. At Quintara, we bring extensive expertise to the field of bispecific antibody development and customization.

## ▶ Service Highlights

### Multiple Types



Asymmetric and symmetric  
Full antibody and antibody fragments

### Large Capacity



1,000 antibody projects/month

### Fast Turnaround




2-3 weeks, from gene to antibody

### Tailored Service





From mL to L  
Optimized plasmid ratio  
Multiple step purification

## ▶ Workflow




# Case Studies




M1: SDS-PAGE Marker  
 Lane 1: Reducing gel (5.00 µg)  
 Lane 2: Non-reducing gel (5.00 µg)  
 Lane 3: BSA (1.00 µg)

Purity (SDS-Page) > 95%

Purity (SEC-HPLC) > 99%




SEC-HPLC




M1: SDS-PAGE Marker  
 Lane 1: BSA (1.50 µg)  
 Lane 2: Non-reducing gel (1.00 µg)  
 Lane 3: Reducing gel (1.00 µg)

Purity (SDS-Page) > 95%

Purity (SEC-HPLC) > 99%




SEC-HPLC



M1: SDS-PAGE Marker  
 Lane 1: BSA (1.00 µg)  
 Lane 2: BsAb (1.60 µg)

Purity (SDS-Page) > 95%

Purity (SEC-HPLC) > 95%



SEC-HPLC